


Digital Image Hub Tutorial

1. Log in with username and password at <http://nbb-slidepath.ncl.ac.uk/dih/login.php>



2. Click on the folder of images

The image shows a screenshot of the 'Digital Image Hub' interface after logging in. The top navigation bar is red and contains links for 'About', 'Contact', 'My Profile', and 'Logout', along with a welcome message 'Welcome HDBR guest (Yuzhu)'. Below the bar, there is a breadcrumb trail 'Location: Browse > Root'. The main content area has a light gray background. It features a 'Browse' tab (which is selected) and a 'Search' button. Below these are two buttons: 'Link' (with a folder icon) and 'Name' (with the value 'H&E and HDBR guest data (Local Image Server)'). A red circle highlights the 'Link' button.

3.



About Contact My Profile Logout

Welcome HDBR guest (Yuzhu)

Location: Browse > Root > H&E and HDBR guest data (Local Image Server)

Link	Name	Stage	Sample number	Slide number	Gene or Stain	Sense? (blank is antisense)	Date Added	Info
<input type="checkbox"/>	Gene Expression Data							
<input type="checkbox"/>	H&E Abnormal Karyotype							
<input checked="" type="checkbox"/>	H&E as a set of representative stages							
<input type="checkbox"/>	H&E in use							

4. Choose the stage of interest



About Contact My Profile Logout

Welcome HDBR guest (Yuzhu)

Location: Browse > Root > H&E and HDBR guest data (Local Image Server) > H&E as a set of representative stages



Link	Name	Stage	Sample number	Slide number	Gene or Stain	Sense? (blank is antisense)	Date Added	Info
<input type="checkbox"/>	CS12							
<input type="checkbox"/>	CS13							
<input type="checkbox"/>	CS14							
<input type="checkbox"/>	CS15							
<input type="checkbox"/>	CS16							
<input checked="" type="checkbox"/>	CS17_N365_Transverse_46XX							
<input type="checkbox"/>	CS18							
<input type="checkbox"/>	CS19_N454_Transverse_46XY							
<input type="checkbox"/>	CS21							
<input type="checkbox"/>	CS23							

5. The folder of images opens (there may be more than one page). Click on the thumbnail to open the image viewer.

[About](#) [Contact](#) [My Profile](#) [Logout](#)

Welcome HDBR guest (Yuzhu)

Location: Browse > Root > H&E and HDBR guest data (Local Image Server) > H&E as a set of representative stages > CS17_N365_Transverse_46XX

<input type="checkbox"/>	Link	Name	Stage	Sample number	Slide number	Gene or Stain	Sense? (blank is antisense)	Date Added	Info
<input type="checkbox"/>		Autotest_0000005833_2013-05-28 17_43_50.scn	CS17	N365	1 2 3 4	H&E		28th May 2013	
<input type="checkbox"/>		Autotest_0000005834_2013-05-28 17_42_39.scn	CS17	N365	5 6 7 8	H&E		28th May 2013	
<input type="checkbox"/>		Autotest_0000005835_2013-05-28 17_42_10.scn	CS17	N365	9 10 11 12	H&E		28th May 2013	
<input type="checkbox"/>		Autotest_0000005836_2013-05-28 17_40_25.scn	CS17	N365	13 14 15 16	H&E		28th May 2013	
<input type="checkbox"/>		Autotest_0000005837_2013-05-28 17_40_08.scn	CS17	N365	17 18 19 20	H&E		28th May 2013	
<input type="checkbox"/>		Autotest_0000005838_2013-05-28 17_37_32.scn	CS17	N365	21 22 23 24	H&E		28th May 2013	
<input type="checkbox"/>		Autotest_0000005839_2013-05-28 17_37_18.scn	CS17	N365	25 26 27 28	H&E		28th May 2013	
<input type="checkbox"/>		Autotest_0000005840_2013-05-28 17_34_37.scn	CS17	N365	29 30 31 32	H&E		28th May 2013	
<input type="checkbox"/>		Autotest_0000005841_2013-05-28 17_34_15.scn	CS17	N365	33 34 35 36	H&E		28th May 2013	
<input type="checkbox"/>		Autotest_0000005842_2013-05-28 17_31_30.scn	CS17	N365	37 38 39 40	H&E		28th May 2013	


 CSV
  PDF
  Images

1
 2
 3
 4
 5
 >
 ||>

6. Use the controls at the bottom of the page to zoom:



Pan across the slide using this tool:




7. Capture a screenshot by clicking the camera icon:



A new window will open. Right click on the image to save it.



8. The thumbnail gallery can be browsed by clicking the arrow icon:



If you have any questions, please contact HDBR@ncl.ac.uk